



Improve your diagnostic testing
in clinical microbiology
with Axonlab

A new benchmark to homogenize
tissue specimens
ProbeAX Evolution



A new benchmark to homogenize tissue specimens for infectious disease diagnostics in disposable tubes

ProbeAX Evolution



Unique and innovative solution

Homogenization through grinding of sample at high speed, releases tissue-associated bacteria, fungi, and biofilms. Suitable for sub-culturing and molecular diagnostics. CE-IVD.



Avoid contamination

Tissue sample is homogenized directly in the sample tube, no additional processing, no contamination risk. Product is Gamma-sterilized.




Ease-of-use

Low noise level during homogenization, added security feature with safety lock during processing.



Saves cost and time

Disposable tubes, fast sample processing, no cleaning of tubes or receptacles.



SpinAX Device

Dispersing and homogenizing system

The SpinAX dispersing and homogenizing system is made to be used in conjunction with the ProbeAX Evolution disposable tubes. The system is designed to grind patient samples and biopsies such as soft tissue, muscle, cartilage- or bone tissue. Aerosol contamination is prevented by the sealed tubes and additional plastic cover over the sample tube during homogenization.

The high-speed grinding system generates an accurate dispersing, grinding and homogenization effect.



The SpinAX is simple to operate with easy options to set time of homogenization. The system has been designed for maximum comfort and runs silently at <56 dB.

ORDER INFORMATION

ARTICLE NUMBER	PRODUCT	PRODUCT DESCRIPTION	PACK CONTENT
AL0509	ProbeAX Evolution	Disposable tube without NaCl, single wrapped	50 pieces
AL0510	ProbeAX Evolution	Disposable tube with 5 ml NaCl (0.85%), single wrapped.	50 pieces
AL0511	ProbeAX Evolution	Disposable tube with 5 ml NaCl (0.85%), double wrapped	50 pieces
AL0511-B	ProbeAX Evolution	Disposable tube with 5 ml NaCl (0.85%), blister double wrapped	50 pieces
AL0500	SpinAX device	Dispersing and homogenizing system, Dimensions (W x D x H): 122 x 178 x 54 mm	1 piece

Visit us – and discover
our products and services
axonlab.com



Headquarters

www.axonlab.com

Axon Lab AG

Täferstrasse 15
CH-5405 Baden-Dättwil
Tel. +41 56 484 80 80
Fax +41 56 484 80 99
info@axonlab.ch

Affiliates

Austria
Belgium
Croatia
Czech Republic
Germany

(for details please refer to www.axonlab.com)

info@axonlab.at
info@axonlab.be
info@axonlab.hr
info@axonlab.cz
info@axonlab.de

Luxembourg info@axonlab.lu
Netherlands info@axonlab.nl
Poland info@axonlab.pl
Slovenia info@axonlab.si
Switzerland (fr) info_f@axonlab.ch

The new generation of Automated Slide Stainer **DAGATRON DX-300**



Designed with safety in mind, it features a specialised ventilation system to protect the technician, while its simple one-button operation and integrated cleaning functions streamline both use and maintenance.

Streamline your process

The DAGATRON Staining System provides a reliable, fully automated solution for Gram and Ziehl/Auramine staining, ensuring consistent results with reduced processing time, reagent usage and labor costs.

Key features

Ventilation System


Ensure the safety and well-being of technicians with the integrated ventilation system, which effectively mitigates potential harm from staining reagent gases.

No Cross Contamination

Slides loaded into separated compartments on the slide tray and stained with new and fresh reagents for each run, eliminating the risk of cross-contamination and ensuring reliable results.

Customization

Select staining protocols based on laboratory needs. Integrate two or three stain protocols into a single system, providing flexibility and adaptability to various laboratory settings.



Enhance efficiency and precision with reliable staining results through the optimized staining protocol.



Stain programs

	Gram	Ziehl-Neelsen
No. of slides	1-20 slides	1-10 slides
Staining time	6 minutes for 10 slides 11 minutes for 20 slides	11 minutes for 10 slides
Reagent	Crystal Violet, Iodine, Acetone Alcohol, Fuchsin and Safranin	Carbol-Fuchsin, Acid Alcohol, Methylene Blue

	Kinyoun	Fluorescence
No. of slides	1-20 slides	1-20 slides
Staining time	14 minutes for 10 slides 19 minutes for 20 slides	10 minutes for 10 slides 15 minutes for 20 slides
Reagent	Carbol-Fuchsin, Acid Alcohol, Methylene Blue	Auramine, Acid Alcohol, Potassium Permanganate

	Wright-Giemsa	May Gruenwald-Giemsa
No. of slides	1-20 slides	1-20 slides
Staining time	6 minutes for 10 slides 9 minutes for 20 slides	15 minutes for 10 slides 17 minutes for 20 slides
Reagent	Wright-Giemsa, Phosphate Buffer	May Gruenwald, Phosphate Buffer, Giemsa

Smart Control System

Customize stain programs to user preferences, with the ability to save up to three programs per stain protocol. Monitor reagent levels and usage in real time for enhanced process control.

Slide Warmer

Improve operational efficiency with the integrated slide warmer. Use for both smear fixation and the Ziehl staining program to improve workflow and reliability.

The new generation of Automated Slide Stainer **DAGATRON DX-300**

Devices

MICROBIOLOGY

ARTICLE NO.	PRODUCT	PACK CONTENTS
DX300-MBSGR00	DX-300 for GRAM	1
DX300-MBSFL00	DX-300 for Fluorescence	1
DX300-MBSZN00	DX-300 for Ziehl-Neelsen	1
DX300-MBSKN00	DX-300 for Kinyoun	1
DX300-MBMAB00	DX-300 for Fluorescence + Ziehl-Neelsen + Kinyoun	1
DX300-MBMGZ00	DX-300 for GRAM + Ziehl-Neelsen + Kinyoun	1
DX300-MBDGF00	DX-300 for GRAM + Fluorescence	1

HEMATOLOGY

ARTICLE NO.	PRODUCT	PACK CONTENTS
DX300-HMSWR00	DX-300 for Wright-Giemsa	1
DX300-HMSMG00	DX-300 for May Grunwald-Giemsa	1

Staining solutions

MYCOBACTERIA

ARTICLE NO.	PRODUCT	PACK CONTENTS
FR-2A(1)	Auramine-Rhodamine	1 x 950 ml
FR-2B(1)	Auramine O	1 x 950 ml
FR-2A(2)	Acid alcohol	1 x 950 ml
FR-2A(3)	Potassium Permanganate	1 x 950 ml
ZR-2A(1)	AFB/Ziehl-Neelsen Carbol fuchsin	1 x 950 ml
ZR-2A(2)	AFB/Ziehl-Neelsen HCL Ethanol	1 x 950 ml
ZR-2A(3)	AFB/Ziehl-Neelsen Methylene blue	1 x 950 ml

GRAM

ARTICLE NO.	PRODUCT	PACK CONTENTS
AL0121_950	AxonGram R1	1 x 950 ml
AL0101_950	AxonGram R2	1 x 950 ml
AL0123_950	AxonGram R3	1 x 950 ml
AL0103_950	AxonGram Fuchsin R4	1 x 950 ml
AL0134_950	AxonGram Safranin R4	1 x 950 ml

HEMATOLOGY

ARTICLE NO.	PRODUCT	PACK CONTENTS
HR-2A(1)	Wright-Giemsa Solution	1 x 950 ml
HR-2A(2)	Phosphate Buffer	1 x 950 ml

Cleaning solution

CHEMICALS

ARTICLE NO.	PRODUCT	PACK CONTENTS
AL0135	Cleaner Ethanol	1 x 950 ml

Controls

MYCOBACTERIA

ARTICLE NO.	PRODUCT	PACK CONTENTS
032-47	Positive	40 slides
032-45	Strongly positive	40 slides
032-40	Simulated sputum	40 slides

GRAM – OTHER

ARTICLE NO.	PRODUCT	PACK CONTENTS
231401	Gram QC	50 slides
032-56	Cryptosporidium spp	20 slides
032-48	Nocardia spp	10 slides



Headquarters

www.axonlab.com

Axon Lab AG

Täferstrasse 15
CH-5405 Baden-Dättwil
Tel. +41 56 484 80 80
Fax +41 56 484 80 99
info@axonlab.ch

Affiliates

Austria info@axonlab.at
Belgium info@axonlab.be
Croatia info@axonlab.hr
Czech Republic info@axonlab.cz
Germany info@axonlab.de

(for details please refer to www.axonlab.com)

Luxembourg info@axonlab.lu
Netherlands info@axonlab.nl
Poland info@axonlab.pl
Slovenia info@axonlab.si
Switzerland (fr) info_f@axonlab.ch

Reagents

 Axon Biotech GmbH
Im Oberfeld 2
DE-94491 Hengersberg
mail@medco.eu, www.medco.eu
A Company of Axonlab

Devices & Reagents

 DAGATRON
117 Tangjung-Ro, Ilsan
Go-Yang, Gyeong-Gi, Korea, 10233
infos@dagatron.com
www.dagatron.com

Automated system for the cultivation of oxygen-sensitive organisms

BACT-R Plus



The BACT-R Plus system replaces all chemical bags for anaerobic, microaerophilic and CO₂ supplementation atmospheres, allowing the best conditions for sustainable use. The major benefits of the BACT-R Plus are its reliability and its capacity to reproduce cultivation conditions over time. All cultivation processes are standardised and complied with the ISO 15189. The traceability system meets the legal requirements of accreditation.

Automated system for the cultivation of oxygen-sensitive organisms

BACT-R Plus

Sustainable

- The reducing valves and the vacuum pump are external and so can be more easily maintained (low maintenance)
- The quality of gas is insured by the gas supplier

Cost-saving

- Jars made in stainless steel adapted to resist to many different pressures
- Three different sizes available for 90/120 mm plates
- Reusable catalysts (till 800 cycles)
- Digital traceability (no paper)



Easy handling

- Direct access to all the programs
- Live follow-up and total control on the screen

Traceability/Quality controls

- All datas are saved on a SD card and can be easily exported
- Barcode and SD readers included
- Quality tests are available (overpressure and underpressure)

ORDER INFORMATION

ARTICLE NO.	PRODUCT	UNIT
7310500TE	BACT-R Plus system (device, vacuum pump and traceability system)	1
7310016E	Ironless Bacto-Jar – 17 plates 90 mm	1
7310014E	Ironless Bacto-Jar – 16 plates 120 mm	1
7310028E	Cristal Bacto-Jar – 17 plates 90 mm	1
7310316AE	Palladium catalyst for anaerobiosis	1



Headquarters

Axon Lab AG

Täferstrasse 15
CH-5405 Baden-Dättwil
Tel. +41 56 484 80 80
Fax +41 56 484 80 99
info@axonlab.ch

Affiliates

Austria
Belgium
Croatia
Czech Republic
Germany

(for details please refer to www.axonlab.com)

info@axonlab.at
info@axonlab.be
info@axonlab.hr
info@axonlab.cz
info@axonlab.de

Luxembourg
Netherlands
Slovenia
Switzerland (fr)

info@axonlab.lu
info@axonlab.nl
info@axonlab.si
info_f@axonlab.ch

www.axonlab.com



SOBIODA
FR-38330 Montbonnot-St-Martin
www.sobioda.com

Optimised Gram stain assay, high stability.

AxonGram



Optimised Gram stain assay, high stability.

AxonGram

Gram's staining of bacteria

Because of the fuchsine the Axonlab's Gram-staining solution employed it stains the bacteria with a higher contrast than other products do. For an optimal colouration, the use of the Axonlab Coloration Sets is recommended. All reagents can of course also be ordered individually.

A clear-cut path towards optimised antibiotics therapy

•Valuable:

A diagnostic aid for every lab

•Simple:


Distinguish bacteria by their
cell wall structure

•Differentiability:

A taxonomic property in
your work

•Fast:

Early identification of
infectious diseases



ORDERING INFORMATION

REAGENTS FOR STAINING SYSTEMS

ARTICLE NO.	PRODUCT	PACK CONTENTS
AL0121	Crystal Violet Solution R1 (0.3% Crystal Violet)	500 ml
AL0101	Lugol's Solution R2	500 ml
AL0123	Decolorizing Solution R3 (94% ethanol; 6% acetone)	500 ml
AL0103	Fuchsine Solution R4	500 ml
AL0134	Safranine Solution R4	500 ml
AL0121_950	Crystal Violet Solution R1 (0.3% Crystal Violet)	950 ml
AL0101_950	Lugol's Solution R2	950 ml
AL0123_950	Decolorizing Solution R3 (94% ethanol; 6% acetone)	950 ml
AL0103_950	Fuchsine Solution R4	950 ml
AL0134_950	Safranine Solution R4	950 ml
AL0111	Object carriers	50 pieces
AL0130	Cleaner Solution	950 ml
AT-3002GZ	Biocentric Borealys	1 device
AL0110V2	Axonlab ColorAX2	1 device

REAGENTS FOR MANUAL USE

ARTICLE NO.	PRODUCT	PACK CONTENTS
AL0100	Crystal Violet Solution R1	500 ml
AL0101	Lugol's Solution R2	500 ml
AL0102	Decolorizing Solution R3	500 ml
AL0103	Fuchsine Solution R4	500 ml
AL0134	Safranine Solution R4	500 ml
AL0111	Object carriers	50 pieces

Axonlab

Headquarters

www.axonlab.com

Axon Lab AG

Täferstrasse 15
CH-5405 Baden-Dättwil
Tel. +41 56 484 80 80
Fax +41 56 484 80 99
info@axonlab.ch

Affiliates

Austria
Belgium
Croatia
Czech Republic
Germany

(for details please refer to www.axonlab.com)

info@axonlab.at
info@axonlab.be
info@axonlab.hr
info@axonlab.cz
info@axonlab.de

Luxembourg info@axonlab.lu
Netherlands info@axonlab.nl
Slovenia info@axonlab.si
Switzerland (fr) info_f@axonlab.ch



Axon Biotech GmbH
Im Oberfeld 2
DE-94491 Hengersberg
mail@medco.eu, www.medco.eu
A Company of Axonlab

Transport and preservation of stool samples – very simple

StorAX stool sample tubes



Transport and preservation of stool samples – very simple

StorAX stool sample tubes

StorAX Native

Transport tube without any medium for the safe transport of stool samples



StorAX Cary-Blair

Transport tubes with Cary-Blair medium, optimised for conserving and detecting intestinal bacteria



StorAX SAF

Transport tubes optimized for preserving parasites in stools and detecting them



The Axonlab range of tubes for stool samples offers with StorAX Native not only the opportunity to transport stools securely and cleanly, but can also be used for preserving and detecting complex pathogenic intestinal bacteria (StorAX Cary-Blair) and intestinal parasites (StorAX SAF). The stool sample tubes, which are CE/IVD marked, are shatterproof (PP), chemically resistant and have a long shelf life.

StorAX Cary-Blair can be stored for 18 months, StorAX SAF for 2 years and StorAX Native even for 4 years.

ORDER INFORMATION

ARTICLE NO.	PRODUCT	PACK CONTENTS
AL0301	StorAX SAF: Transport medium for conserving parasites (protozoa and worm eggs) in stools	25 stool tubes with SAF solution
AL0301S	StorAX SAF Single: Transport medium for conserving parasites (protozoa and worm eggs) in stools	100 stool tubes with SAF solution packaged individually
AL0302	StorAX Cary-Blair: Transport medium for conserving intestinal bacteria	25 stool tubes with a Cary-Blair solution
AL0302S	StorAX Cary-Blair Single: Transport medium for conserving intestinal bacteria	100 stool tubes with Cary-Blair solution packaged individually
AL0303	StorAX Nativ: Transport medium for raw stool sample	25 stool tubes without medium
AL0303S	StorAX Nativ Single: Transport medium for raw stool sample	100 stool tubes without medium packaged individually
AL0310	Filtration unit	25 filtration units and 25 centrifuge tubes
AL0311	StorAX SET SAF with filtration unit and centrifuge	25 stool tubes StorAX SAF 7.5 ml 25 filtration funnels and 25 centrifuge tubes

StorAX Nativ – Detection

Adenoviruses
Calprotectin
Entamoeba histolytica antigen
Helicobacter pylori antigen
Leukocytes
Mycobacteria
Viruses
Occult blood
Rotaviruses, Noroviruses, Adenoviruses
Mold

StorAX Cary-Blair – Detection

Stool bacteriology, Salmonella, Shigella, Campylobacter
Yersinia spp. / *Aeromonas spp.*
Clostridium difficile
Enterotoxigenic *E. coli* (ETEC)
Verocytotoxin-producing *E. coli* (EHEC/STEC)
Vibrio cholerae

StorAX SAF – Detection

Amoebae
Blastocystis hominis
Giardia lamblia antigen
Cryptosporidia
Protozoa
Strongyloides
Worm eggs



Headquarters

Axon Lab AG
Täferstrasse 15
CH-5405 Baden-Dättwil
Tel. +41 56 484 80 80
Fax +41 56 484 80 99
www.axonlab.com
info@axonlab.ch

Affiliates

Austria info@axonlab.at
Belgium info@axonlab.be
Croatia info@axonlab.hr
Czech Republic info@axonlab.cz
Germany info@axonlab.de

(for details please refer to www.axonlab.com)

Luxembourg info@axonlab.lu
Netherlands info@axonlab.nl
Slovenia info@axonlab.si
Switzerland (fr) info_f@axonlab.ch



Agar plates for the Kirby-Bauer
sensitivity testing –
square (120 mm) and round (ø 90 mm)
Mueller-Hinton agar plates




Agar plates for the Kirby-Bauer sensitivity testing – square (120 mm) and round (ø 90 mm)

Mueller-Hinton agar plates

Mueller-Hinton agar is a standard medium used for antimicrobial sensitivity testing using the agar diffusion test and also as a complete medium for the cultivation of gram-positive and/or gram-negative bacteria. The composition offers favourable conditions for growth on the one hand, and the widespread absence of sulphonamide antagonists on the other hand.

The round and square plates are suitable for manual methods and for automated readers and incubators, such as SIRscan Orion® or SIRscan LYNX® from i2a.

Square plates (120 mm) make it possible to apply up to 16 antibiotic discs and thus reduce the number of agar plates for antibiograms and imaging documentation.



ORDERING INFORMATION

ARTICLE NO.	PRODUCT	STORAGE	PACKAGE CONTENTS
AL0901	Mueller-Hinton agar plates (120 mm)	+4 to +8°C	50 square Petri dishes
AL0909	Mueller-Hinton agar plates (90 mm)	+4 to +8°C	20 round Petri dishes

PRODUCT ACCESSORIES SPECIFICATIONS

ARTICLE NO.	PRODUCT DESCRIPTION	STORAGE	PACKAGE CONTENTS
AL0604	NaCl 0.85% solution, CE-IVD	15°C to 30°C	100 x 3 ml
On demand	ATCC control strains		
On demand	Antibiotic disks from i2a		
800256	Square dispenser from i2a		
800226	Round dispenser from i2a		



Headquarters

Axon Lab AG

Täferstrasse 15
CH-5405 Baden-Dättwil
Tel. +41 56 484 80 80
Fax +41 56 484 80 99
info@axonlab.ch

www.axonlab.com

Affiliates

Austria
Belgium
Croatia
Czech Republic
Germany

(for details please refer to www.axonlab.com)

info@axonlab.at
info@axonlab.be
info@axonlab.hr
info@axonlab.cz
info@axonlab.de

Luxembourg
Netherlands
Slovenia
Switzerland (fr)

info@axonlab.lu
info@axonlab.nl
info@axonlab.si
info_f@axonlab.ch



Axon Biotech GmbH
Im Oberfeld 2
DE-94491 Hengersberg
mail@medco.eu, www.medco.eu
A Company of Axonlab

Ready-to-use solution
for your identifications and
antibiotic sensitivity tests
NaCl tubes and DensitomAX



Ready-to-use solution for your identifications and antibiotic sensitivity tests

NaCl tubes and DensitomAX



Standardise your suspension
More efficiency



Save time
Increase productivity



No contaminations
Reliable results



McFarland suspension
preparation
Coordinated total solution

NaCl tubes 0.45% 3 ml **Ready-to-use solutions**

Ready-to-use solution for pathogen
identification and/or resistance testing.


These tubes are suitable for commercially
available devices.



NaCl tubes 0.85% 2 ml **Ready for manual AST**

Ready-to-use solution for resistance testing
by agar diffusion test.

These tubes are suitable for manual
AST methods.



With the included adapters, tubes with diameters of 12, 16 and 18 mm can be used.



DensitomAX

No additional densitometer needed in your lab

DensitomAX offers you a total solution for your McFarland suspension preparation. DensitomAX is a compact and efficient bench top densitometer, factory calibrated with our own NaCl test tubes to ensure high accuracy of the results. DensitomAX is designed to measure turbidity in the range of 0.00 to 6.0 McFarland

units. If your laboratory uses different type of tubes, DensitomAX can handle also other tube types (12 mm, 16 mm and 18 mm diameter, round or flat bottom, glass or plastic). DensitomAX can also be used in other life science applications where turbidity measurements are needed.

ORDER INFORMATION

ARTICLE NUMBER	PRODUCT DESCRIPTION	PACK CONTENT
AL0603	NaCl 0.45% solution / Storage: 15°C to 30°C / CE-IVD	100 x 3 ml
AL0604	NaCl 0.85% solution / Storage: 15°C to 30°C / CE-IVD	100 x 2 ml
BS-050104-AXW	DensitomAX / CE	1 unit

**Do you need a different tube format or filling volume for your workflow?
Please contact us.**

Visit us – and discover
our products and services
www.axonlab.com



Headquarters

Axon Lab AG

Täferenstrasse 15
CH-5405 Baden-Dättwil
Tel. +41 56 484 80 80
Fax +41 56 484 80 99
info@axonlab.ch

www.axonlab.com

Affiliates

Austria
Belgium
Croatia
Czech Republic
Germany

(for details please refer to www.axonlab.com)

info@axonlab.at
info@axonlab.be
info@axonlab.hr
info@axonlab.cz
info@axonlab.de

Luxembourg info@axonlab.lu
Netherlands info@axonlab.nl
Slovenia info@axonlab.si
Switzerland (fr) info_f@axonlab.ch

NaCl tubes



Axon Biotech GmbH
Im Oberfeld 2
DE-94491 Hengersberg
mail@medco.eu, www.medco.eu
A Company of Axonlab



DensitomAX
SIA Biosan
LV-1067, Riga, Latvia
www.biosan.lv



UrinAX from Axonlab

Bacteriology of urinary tract
infections in one simple step

UrinAX from Axonlab

Bacteriology of urinary tract infections in one simple step

Axonlab is setting new standards for the bacteriology of urinary tract infections. The research and development department at Axonlab, the manufacturer of the UrinAX product range, has developed five different immersion media, tailored to meet individual customers' needs.

These not only allow for bacterial counts and the selective growth of bacteria in urine (MacConkey, Enterococcus, Cetrимide and Malt agar), but also give an indication after just 18-24 hours if Trimethoprim/Sulfamethoxazole resistant bacteria are present in the urine. This makes it possible to target antibiotic treatment more accurately.


Urine dip slides, transport medium, bacterial counts and selective growth of bacteria in urine

- Excellent bacteria growth, including Streptococcus B
- Extra-large CLED agar surface makes assessing the results and transferring of colonies easier
- MacConkey selective agar for gram-negative bacteria
- Selective agar for *Pseudomonas aeruginosa* (Cetrимide agar), malt agar for yeast and selective agar for detecting enterococci

UrinAX TS:

From the bacteria count to the selective growth of bacteria and an indication of the antibiotic resistance in the urine takes only 24 hours.

	Conventional bacteriology	UrinAX TS
Bacterial count	24h	24h
Resistance	24h	
Total time until result	48h	24h



Mini Incubator

Compact size, maximum performance.

The Mini Incubator is the ideal solution for labs. The interior and the housing are robust and corrosion-resistant. Because of its small size, the Mini Incubator will fit in any laboratory.



Stand

The Plexiglas stand is designed for daily use in the lab. It holds 10 UrinAX dip slides.



PRODUCT SPECIFICATIONS

PRODUCT	AGAR TYPES	PACK CONTENTS	STORAGE	ARTICLE NO.
UrinAX CL/MC/TS	CLED/MacConkey/Trimethoprim-Sulfamethoxazole (TS)	10	15-25°C	AL0201
UrinAX CL/MC	CLED/MacConkey	10	15-25°C	AL0202
UrinAX CL/MC/E	CLED/MacConkey/Enterococcus	10	15-25°C	AL0203
UrinAX CL/MC/CETRI	CLED/MacConkey/Cetrimide	10	15-25°C	AL0205
UrinAX CL/MC/MALT	CLED/MacConkey/Malt	10	15-25°C	AL0206

SUPPLEMENTARY PRODUCTS

PRODUCT	PRODUCT DESCRIPTION	PACK CONTENTS	ARTICLE NO.
Mini Incubator	Dimensions (W x L x H): 28.5 x 28 x 33.5 cm, 230 V	1 device	AL5110
Stand	Plexiglas, capacity: 10 UrinAX dip slides	1 item	AL02 Stand

UrinAX from Axonlab

Bacteriology of urinary tract infections in one simple step



- The broad range of UrinAX products provides a simple method for carrying out bacterial counts and the selective growth of bacteria in urine
- The large CLED agar surface makes assessing the results and transferring colonies easier
- Infection analysis within 18-24 hours



Headquarters

www.axonlab.com

Axon Lab AG

Täferenstrasse 15
CH-5405 Baden-Dättwil
Tel. +41 56 484 80 80
Fax +41 56 484 80 99
info@axonlab.ch

Affiliates

Austria
Belgium
Croatia
Czech Republic
Germany

(for details please refer to www.axonlab.com)

info@axonlab.at
info@axonlab.be
info@axonlab.hr
info@axonlab.cz
info@axonlab.de

Luxembourg info@axonlab.lu
Netherlands info@axonlab.nl
Slovenia info@axonlab.si
Switzerland (fr) info_f@axonlab.ch



Axon Biotech GmbH
Im Oberfeld 2
DE-94491 Hengersberg
mail@medco.eu, www.medco.eu
A Company of Axonlab

Vacuum solutions for all
applications necessitating
the safe disposal of liquids
Mini Vacuum Pump



A safe disposal of biological liquids


Mini Vacuum Pump

Mini-Vac power combines all advantages of aspiration, disposal and filtration of biological and chemical fluids. The Mini Vacuum Pump features a powerful vacuum, space saving footprint and a low noise level. Aspiration of fluids, such as cell culture media, cell suspensions, supernatants etc., disposal of biological and chemical waste fluids.

- Very strong pump (9.5 l / min (current pump ± 7 l / min))
- Vacuum filtration and degassing of fluids
- Powerful vacuum, adjustable from 0 to -900 mbar
- Very effective noise insulation (< 60 dB)
- All parts in contact with fluids are autoclavable
- Adjustable vacuum for gentle aspiration of small volumes
- Innovative design
- Easy to clean and to carry (handle)
- Possibility of installing a second Vacuum-Hand handle onto the same device
- External level sensor prevents damage to the pump and abolishes the need to clean the sensor.

HAS-1 Hand operator set

Axonlab also offers practical HAS-1 Hand operator sets for easy handling and convenient aspirating.



- 1 Single channel adapter for 200 μ l tips
- 2 8-channel adapter for 200 μ tips
- 3 8-channel adapter with 35 mm stainless steel pin
- 4 Single channel adapter with 40 mm stainless steel pin
- 5 Single channel adapter with 125 mm stainless steel pin
- 6 SH-6, Stand for HAS-1
- 7 Handheld vacuum controller

ORDER INFORMATION

ART.-NO.	REF.-NO.	PRODUCT	PACK CONTENTS
11005031	M30276	Mini Vacuum system, instruction manual	1 device
11037004	M30276R	Mini Vacuum system with accessories. Level sensor, vacuum bottle (PP, 2l), filling/venting closure, tubing set, PTFE filter, instruction manual	1 complete set
12055681	BS040118PK	HAS-1 Hand operator set with Stand SH-6, 8-channel adapter, single channel adapter and handheld vacuum controller	1 kit
11005741	VF50ANPPT004AC01	Vent filter Ø 50 mm barb connectors PP transparent membr. PTFE 0.45 μ m	5 pieces
11037005	M30276S	Level Sensor	1 device
12045484	M30378	Foot Switch for Mini-Vac	1 device

PRODUCT SPECIFICATION

MINI-VAC POWER PUMP	
Voltage	230 V/50 Hz
Output	60 W
Working tube	Silicone tube, 0 6/10 mm, length 1'000 mm
Vacuum/aspiration power	0 to -900 mbar
Vacuum bottle	Nalgene heavy duty vacuum bottle PP (2 L), autoclavable
Closure	Nalgene filling/venting cap with two fittings, autoclavable
Safety filter	1 Vent filter Ø 50 mm barb connectors PP transparent membr. PTFE 0.45 μ m GVS
Net Weight	4.5 kg



Headquarters

Axon Lab AG
Täferstrasse 15
CH-5405 Baden-Dättwil
Tel. +41 56 484 80 80
Fax +41 56 484 80 99
www.axonlab.com info@axonlab.ch

Affiliates

Austria info@axonlab.at
Belgium info@axonlab.be
Croatia info@axonlab.hr
Czech Republic info@axonlab.cz
Germany info@axonlab.de

(for details please refer to www.axonlab.com)

Luxembourg info@axonlab.lu
Netherlands info@axonlab.nl
Slovenia info@axonlab.si
Switzerland (fr) info_f@axonlab.ch

Enrichment medium
for EHEC bacteria
for clinical microbiology
EHEC enrichment broth



Enrichment medium for EHEC bacteria for clinical microbiology **EHEC enrichment broth**

**The EHEC enrichment broth is an enrichment medium
for EHEC bacteria directly from stool samples.
Also detects weak toxin creators.**

A fundamental problem when isolating EHEC bacteria from stool samples is that apathogenic *E. coli*, which are resident in the gut, are always present at the same time as EHEC bacteria but in larger numbers and with a wide range of types. Enrichment in a suitable culture medium is therefore necessary to isolate suspected EHEC bacteria. Yet a pre-enrichment or enrichment of the relevant bacteria has also been found to be necessary for finding evidence of Shiga toxins.

Bacterial growth and the formation of Shiga toxins by EHEC bacteria are stimulated in the EHEC enrichment broth. After the enrichment Shiga toxins can be detected using an ELISA.



REAGENTS

COMPOSITION PER LITRE	
Peptone	20 g
Sodium chloride	5 g
Antibiotic mix	approx. 16 mg
Growth factors	200 mg

100 test tubes with 4.5 ml of EHEC enrichment broth are delivered with each pack.

ORDER INFORMATION

ARTICLE NO.	PRODUCT	STORAGE	PACK CONTENTS
AL0602	EHEC enrichment broth: Enrichment of EHEC from stool samples	2°C to 8°C	100 x 4.5 ml

SUPPLEMENTARY PRODUCTS

ARTICLE NO.	PRODUCT	STORAGE	PACK CONTENTS
AL0302	StorAX Cary-Blair: Transport system for the preservation of intestinal bacteria	15°C to 30°C	25 stool test tubes with Cary-Blair solution
AL0303	StorAX Nativ: Transport system for stool	15°C to 30°C	25 stool test tubes without medium



Headquarters

Axon Lab AG
Täferstrasse 15
CH-5405 Baden-Dättwil
Tel. +41 56 484 80 80
Fax +41 56 484 80 99
www.axonlab.com
info@axonlab.ch

Affiliates

Austria info@axonlab.at
Belgium info@axonlab.be
Croatia info@axonlab.hr
Czech Republic info@axonlab.cz
Germany info@axonlab.de

(for details please refer to www.axonlab.com)

Luxembourg info@axonlab.lu
Netherlands info@axonlab.nl
Slovenia info@axonlab.si
Switzerland (fr) info_f@axonlab.ch



Axon Biotech GmbH
Im Oberfeld 2
DE-94491 Hengersberg
mail@medco.eu, www.medco.eu
A Company of Axonlab



Safety SubCulture Unit 2

Versatile dispensing unit for a safer subculture technique

Make your lab safer and more efficient with Safety SubCulture Unit 2.

It enables simple and safe needleless transfer of samples from culture bottles to slides and media plates.

Designed to fit securely onto a range of culture bottles, the plastic dispensing tip makes it easy to dispense samples onto slides and plates.

The protective cap captures gas and culture discharge, reducing potential exposure.


The new standard for safety in blood culture sampling



Choose the safer solution for subculture

Safety SubCulture Unit 2:

- Needleless dispensing of blood culture samples, reducing the risk of sharps injury
- Secure fit onto culture bottle septum
- Protective cap over the dispensing tip reduces potential exposure
- Safe and easy sample dispensing procedure onto slides and plates
- Fits a wide variety of culture bottles



Visit itlbiomedical.com to request your sample today

Product Code:

A100722: Sterile, 288/case, 48/box, 6 boxes/case

Product Information:



Not Made with Natural Rubber Latex



Single Use Only





Safety SubCulture Unit

Needleless, efficient sample transfer from blood culture bottles


Make your lab safer and more efficient with the Safety SubCulture Unit. It enables safe, needleless transfer of samples from positive blood culture bottles.

This IVD blood culture device has been developed to safely transfer positive samples using the drop method or a needleless syringe to transfer measured samples.

Designed to fit with a range of blood culture bottles, the plastic safety tip pierces the bottle septum and makes it easy to insert and remove.

The integrated filter cap also captures gas and culture discharge, reducing potential exposure.

The new standard for safety in blood culture sampling



Choose the safer solution for subculture

Safety SubCulture Unit:

- Two controlled sampling options
 - Dispense drops directly from culture bottle
 - Needleless syringe withdrawal and dispensing
- Integrated filter cap to capture gas and culture discharge, reducing exposure
- Plastic safety tip pierces the bottle septum and makes it easy to insert and remove from the culture bottle
- Fits a wide variety of culture bottles



Two controlled sampling options:

1. Built-in sampling channel to dispense drops directly from culture bottle.
2. Female luer connection for needleless syringe withdrawal and dispensing.

1



2



Visit itlbiomedical.com to request your sample today

Product Code:

A100720: Sterile, 360/case, 60/box, 6 boxes/case

Product Information:



Not Made with Natural Rubber Latex



Single Use Only





SampLok® Adapter Cap 3

Small, simple and efficient blood culture sample collection

Designed to create minimal biohazard sharps waste, the compact SampLok® Adapter Cap 3 is a versatile device for safe blood sample collection into culture bottles and vacuum test tubes.

Versatile and compact, SampLok® Adapter Cap 3 fits securely on a range of culture bottles and vacuum tubes and ensures the sampling needle is well positioned in the culture bottle septum or tube cap.

The smaller size of these adapter caps reduces biohazard sharps waste, enhancing the efficiency and sustainability of your blood culture sample collection process.

The new standard for safety in blood culture sampling



Improve blood culture sample collection safety, productivity and efficiency

SampLok® Adapter Cap 3:

- Compact size generates reduced biohazard sharps waste
- Secure fit makes it easy to handle while inserting, filling and removing culture bottles and vacuum tubes
- Flexible design fits with a variety of culture bottles and vacuum test tubes
- Optional insert for use with sample tubes and long-neck culture bottles
- Can be attached to blood collection and winged needle sets



Visit itlbiomedical.com to request your sample today

Product Code:

A100760: Non-sterile, 600/case, 100/box, 6 boxes/case

A100761: Insert, non-sterile, 840/case, 70/box, 12 boxes/case

Other configurations available, visit our website for more information

Product Information:



Not Made with Natural Rubber Latex



Single Use Only



1.888.411.2851 sales@itlbiomedical.com itlbiomedical.com

SampLok® is a registered trademark of ITL Corporation, Canberra, Australia.
Copyright © 2018 ITL Corporation. All rights reserved. ITL-BR-0067 Rev01